For veterinary use only

INTENDED USE

The Vcheck SDMA is an *in vitro* diagnostic test kit for the quantitative measurement of SDMA concentration in canine or feline serum and plasma. Symmetric dimethylarginine (SDMA) is a methylated derivative of arginine produced from intranuclear methylation of L-arginine residues. SDMA is primarily excreted by the kidneys, making it an endogenous marker of renal function. SDMA is a more reliable indicator of kidney function than creatinine and the measurement of SDMA can be used to detect both acute kidney injury (AKI) and chronic kidney disease (CKD) earlier than ever. The BIONOTE Vcheck SDMA is designed to be used only by veterinarians.

PRINCIPLE

This product uses an anti-SDMA antibody that specifically binds to SDMA, using fluorescence measurement and immunoassay for SDMA measurement in serum or plasma (heparin). After mixing the sample with Pretreatment Buffer, the supernatant is recovered by centrifugation, diluted with assay diluent, and mixed with a tablet made of anti-SDMA antibody bound to fluorescent particles. Then, when this sample mixture is dropped on the sample hole of the SDMA test device, the sample mixture moves along the membrane inside the device along with the capillary phenomenon, and is adsorbed onto the membrane area coated with SDMA. As a result, when SDMA activates the fluorescent material after binding to ultraviolet rays through the competitive reaction with SDMA coated on the membrane, the reflected light is read to measure the SDMA concentration in serum or plasma. Immunofluorescent assay is performed using the optical system inside the BIONOTE Vcheck analyzer, and the SDMA concentration in the sample is displayed by an algorithm for the SDMA result value.

MATERIALS PROVIDED

Reagent	5 Tests/Kit	10 Tests/Kit
① Vcheck SDMA Test device	5	10
② Assay diluent bottle	1	1
③ Pretreatment Buffer	2	2
④ Disposable tablet pipette (orange bulb)	5	10
⑤ Disposable pipette tip	20	40
⑥ 1.5 mL Tube	10	20
⑦ Instructions for use	1	1

MATERIALS REQUIRED, BUT NOT PROVIDED

1. BIONOTE Vcheck Analyzer
2. 25 µl pipette
3. 50 µl pipette
4. 100 µl pipette

STORAGE AND STABILITY

1. Store the test kit at 2-8 °C. **DO NOT FREEZE.**
2. Do not store the test kit in the direct sunlight.
3. The test kit is stable until the expiry date that is marked on the package label.

Reagent	Open status	Storage	Stability	Note
Test device and tablet pipette	Unopened	2-8 °C, Sealed	9 months	Finished product
	Opened	Do not store	-	Use directly
Assay diluent	Unopened	2-8 °C, Sealed	9 months	Finished product
	Opened	2-8 °C, Resealed	9 months	-
Pretreatment Buffer	Unopened	2-8 °C, Sealed	9 months	Finished product
	Opened	2-8 °C, Resealed	9 months	-

PRECAUTIONS

1. The test kit is for canine or feline use only. Do not use for other animals.
2. The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the aluminum foil pouch.
3. Do not reuse test components.
4. Do not touch the membrane in the result window of the test device.
5. Do not use the test kit beyond the stated expiry date marked on the label.
6. Strictly follow the test procedure (E.g. adequate sample volume), as failure to do so may adversely affect test performance and/or produce invalid results.
7. Although the Vcheck SDMA Test Kit offers simple and quick quantitative measurement of SDMA concentration in canine or feline serum and plasma, there may be a difference in the detection performance with other clinical or laboratory methods with more sophisticated principles.
8. Do not use the test kit if the pouch is damaged or the seal is broken.
9. Do not mix components from different lot numbers, the components in this kit have been quality control tested as a standard batch unit.

10. All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
11. Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulations.
12. Hemolyzed samples or serum with precipitate such as fibrin may give erroneous results.
13. Pretreatment Buffer is acidic, so care should be taken not to touch the skin.
14. Available only in countries other than the United States and Europe.
15. Professional veterinarian should make a final diagnosis based on the results of this product, other test results and clinical findings.
16. BIONOTE Vcheck analyzer is recommended to use at 15-30 °C.

COLLECTION AND PREPARATION OF SAMPLE


1. Serum or plasma of both feline and canine should be used with this test. A method for preparing the specimen is as follows.
 2. **[Serum]** Collect the whole blood into a blood collection tube containing **NO anticoagulant**. Leave to settle for 30 minutes for blood coagulation and then centrifuge to obtain the serum supernatant.
 3. **[Plasma]** Collect the whole blood into the collection tube containing anticoagulant (only heparin) and then centrifuge to obtain plasma.
- If serum and plasma samples are not tested immediately, they should be refrigerated at 2-8 °C and used within 2 weeks. For longer storage, serum or plasma can be frozen (-20 °C or colder).

TEST PROCEDURE

- * Allow all kit components and samples to reach room temperature (15-30 °C) for at least 30 minutes prior to testing.
- * Prepare the same number of the test device, anticoagulants, assay diluent and 1.5 mL tube as the number of the samples.

[Coding]

1. Turn on V200 Analyzer and select "Standard Test".


2. Remove the test device from the aluminum foil pouch. Once the "Insert Device" is displayed in the screen, insert the test device.


[Pretreatment]

3. Use a 100 µl pipette to draw 100 µl of the sample and add to a 1.5 mL tube.


4. Use a 25 µl pipette to draw 25 µl of the Pretreatment Buffer and add it to the 1.5 mL tube containing the sample.

- * The pipette tip will be clogged as soon as the Pretreatment Buffer touches the sample. The end of pipette tip should be applied onto the inner wall of the tube.
- * Be careful not to touch the skin with the Pretreatment Buffer.


5. Tap with fingers 6-8 times or vortex to mix well.

- * Ensure that the sample and Pretreatment Buffer are well mixed to form white sediment.


6. Centrifuge the well-mixed solution at 6,000 x g for 5 minutes.

- * At 10,000 RPM for 5 minutes with the centrifuge provided by BioNote.


[Dilution of sample & Measurement]

7. Collect 50 µl of the separated supernatant and transfer to a NEW 1.5 mL tube.

- * Be careful not to bring the sediments at the bottom when collecting the supernatant because the sediment is not fixed.

8. Use 50 µl pipette to collect 50 µl of the assay diluent and add it to the tube containing the supernatant.

- * Be sure to mix the supernatant and the diluent thoroughly.

9. Use a disposable tablet pipette to mix well by squeezing the rubber bulb gently several times until the white tablet dissolves in the pipette completely.

- * Squeeze the orange bulb gently to minimize foam.

10. Add all of the mixed sample/tablet into the sample hole of the test device using the disposable tablet pipette and press the [START] to initiate testing.

11. The V200 analyzer will display the test result on the screen after 5 minutes.

12. Remove the test device.

- * Strictly follow the test procedure including the amount of sample (100 µl) used and the test time (5 min), as failure to do so may adversely affect test performance and/or produce invalid results.

INTERPRETATION OF THE RESULT

1. Read the concentration value of SDMA appearing on the display of the BIONOTE Vcheck Analyzer. (10-100 µg/dL)
2. If "↓ 10 µg/dL" appears on the display, it means the concentration of SDMA the specimen is less than 10 µg/dL.
3. If "↑ 100 µg/dL" appears on the display, it means the concentration of SDMA the specimen is greater than 100 µg/dL.
4. If the [Invalid] result appears on the screen, a retest shall be carried out.

REFERENCE RANGE

≤ 14 µg/dL	14.1 - 19.9 µg/dL	≥ 20 µg/dL
Normal (≤ 16 µg/dL in puppies*)	Suspected (Check other evidence of kidney disease. If there are no other symptoms, a retest is recommended in 2-3 weeks.)	Kidney disease probable

* Mildly increased SDMA concentrations (14-16 µg/dL) in puppies should be interpreted in light of the growth phase as well as other evidence of kidney disease.

PERFORMANCE CHARACTERISTIC

1. Measuring Range
10-100 µg/dL of canine or feline SDMA is measurable in serum or plasma. For ↑ 100 µg/dL results, samples may be diluted with negative sample and measured again for a numerical result. To calculate the final canine or feline SDMA concentration, the measured result must be multiplied by the dilution factor.

SCREEN MESSAGES AND TROUBLE SHOOTING

[V200]

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. Solution: Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. Solution: Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. Solution: Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. Solution: Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed. Solution: Reconnect the analyzer and external device.
Barcode Error	If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Extremely High Total Hemoglobin	The measured total hemoglobin is out of the range of 7 to 23 g/dL. Solution: This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	The test is invalid. Solution: Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Calibration Overdue	The calibration is overdue. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Not Supported Device	A test device that is not supported by the analyzer has been loaded. Solution: Check whether the test device is manufactured by BioNote, Inc.
EEE	Internal error has occurred. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.

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