

Caniv-4 (Babesia)

General Information

Purpose : Detection of *D. immitis* antigen, canine antibody against *E. canis*, *B. gibsoni*, and *A. phagocytophilum/A. platys*

Principle : Immunochromatographic assay


Specimen : Canine whole blood, plasma or serum

Materials provided :


- Anigen Rapid Caniv-4 (Babesia) test device
- Assay diluents bottle
- Anticoagulant tube
- Disposable capillary tube

Special Features :

- Simultaneous detection of four kinds of Vector Borne diseases
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity



Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RC2133DD	Caniv-4 (Babesia)	Quattro Device	10 Tests/Kit