

Anigen Rapid FPV Ag

General Information

Purpose : Detection of Feline Parvovirus antigen

Principle : Immunochromatographic assay

Specimen : Feline feces

Materials provided :

- Anigen Rapid FPV test device
- Assay diluents tube
- Disposable swab
- Disposable dropper




Special Features :

- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity
 - Sensitivity : 97% vs. PCR
 - Specificity : 98.5% vs. PCR



Test Procedure

[Proper sample amount collected by a swab]



Ordering Information

Cat. No.	Description	Type	Pack size
RG1203DD	FPV Ag	Device	10 Tests/Kit